



## Controls, lyophilised

### Instructions for use

The lyophilisates must be re-suspended before use, serve as control substances for the inhibitor analysis and are used e.g. in the Brilliant Black Reduction Test (BRT).

### Composition and concentration

The lyophilisates contain after re-suspending the substance-concentration mentioned on the vial +/- 10 %.

### Storage and use by date

The lyophilisates can be conserved at 6 – 10 °C, dark, until the mentioned date on the packing.

### Instructions for use



- 1) Flip and tear off aluminium seal of the vial.
- 2) Remove the buckler from the vial and lay it up side down on the bench.
- 3) Pipette 4.5. ml resp. 15 ml deionized, sterile (20 – 30 °C) into the vial (recommendation art.-no. 9202 res p. 9203, water deionized sterile, 4.5 ml resp. 15 ml).
- 4) Close vial with the buckler.
- 5) Shake the vial or leave it on the bench. As soon as the lyophilisate is re-suspended completely, the solution is ready for use.
- 6) Store the solution refrigerated (6 – 10 °C) and use it max. three days.
- 7) In case of contamination the solution may be rendered unusable or may lead to incorrect results.
- 8) The solution may be portioned in small volumes for daily demand and frozen (– 20 °C) for later use, labelling is recommended.

### Safety instructions

Hazard labelling is not necessary according to EC Directive 67/648/EEC. Read material safety data sheet.

### Certificate of analysis

Every lot is checked internally and externally, certificates of analysis are available.